



The **Brightest Ideas** In Dental Implants





Dear Customers and Colleagues

Today, dental implants have become an indispensable part of dental treatment options. With the globalization of medical infrastructures and higher standards of living, implant applications have rapidly become common.

Southern Implants has been a manufacturer and distributor of dental implants since 1987. Today, the Southern group is recognized as a leading bio-medical engineering entity, with major intellectual property and capabilities in implantable devices, arthroplasties, tissue regeneration, stem cells and cryoscience. The top-end professional users, who want more choices, have driven the product range expansion to enormous and exciting heights. Striving for excellence and meeting customer needs has lead to our wide product range characterized by numerous unique and innovative products which include:

- 3 interfaces: External Hex, Internal morse taper/octagon, and Tri-Nex.
- Many products optimized for primary stability and suited for immediate loading.
- The only angled-top tapered screw-form 12°, 24° and 36° Co-Axis implant.
- Implant lengths from 6mm to 20mm and diameters from 2.90mm to 10mm.
- A surface which continues to out-perform that which it is trailed against.
- Color-coded components for easy part recognition.
- 55° Zygomatic implant, optimized for load distribution.
- Compatibility with major brands, giving the patient more options.
- The MAX, wide diameter implant for molar teeth replacement.

Striving for excellence is synonymous with the search to improve. At Southern the development starts with computer simulation and finite element modeling. This is followed by extensive laboratory trials and testing. Finally, clinical research has taken on a new dimension in our overall strategy where our preference is for independent RCTs.

Our sincere thanks to all specialists, dentists and technicians who give continual feedback, suggestions and input. The products here are our interpretation of your needs.

Yours sincerely



Graham Blackbeard Managing Director Southern Implants





















Contents

Intro	oduction and Welcoming Letter	Inside Front Cover	
	Ø3.25mm IBN & IBNT Implants, Components and Drilling Information	Page 04 Page 05	
	Ø3.30mm INPi Implant, Components and Drilling Information INPi Prosthetics Flowchart & Additional Information	Page 06 Page 07	
	Ø3.75mm & Ø4.00mm IB & IBT and IBPS Implants, Components and Drilling Information. IB, IBT & IBPS Prosthetics Flowchart.	Page 08 Page 09	
	Ø5.00mm BA & BAT Implants, Components and Drilling Information. BA & BAT Prosthetics Flowchart.	Page 10 Page 11	
	Ø6.00mm BBB & BBBT Implants, Components and Drilling Information	Page 12 Page 13	
MA	MAX Implant Introduction		
	MAX-7 Implant, Components and Drilling Information. MAX-7 Prosthetic Flowchart.	Page 16 Page 17	
	MAX-8 Implant, Components and Drilling Information. MAX-8 Prosthetic Flowchart.	Page 18 Page 19	
	MAX-9 Implant, Components and Drilling Information	Page 20 Page 21	
Co-	Axis Implant Introduction	Page 22-23	
Co-	Axis Implant Introduction. Ø4.00mm IBT12d Co-Axis Implant, Components and Drilling Information. IBT12d Prosthetic Flowchart.	Page 22-23 Page 24 Page 25	
Co-	Ø4.00mm IBT12d Co-Axis Implant, Components and Drilling Information.	Page 24	
	Ø4.00mm IBT12d Co-Axis Implant, Components and Drilling Information. IBT12d Prosthetic Flowchart. Ø5.00mm BAT12d & BAT24d Co-Axis Implants, Components and Drilling Information.	Page 24 Page 25 Page 26	
	 Ø4.00mm IBT12d Co-Axis Implant, Components and Drilling Information. IBT12d Prosthetic Flowchart. Ø5.00mm BAT12d & BAT24d Co-Axis Implants, Components and Drilling Information. BAT12d and BAT24d Prosthetic Flowchart. Ø5.00mm BAT36d Co-Axis Implants, Components and Drilling Information. 	Page 24 Page 25 Page 26 Page 27 Page 28	
	 Ø4.00mm IBT12d Co-Axis Implant, Components and Drilling Information. IBT12d Prosthetic Flowchart. Ø5.00mm BAT12d & BAT24d Co-Axis Implants, Components and Drilling Information. BAT12d and BAT24d Prosthetic Flowchart. Ø5.00mm BAT36d Co-Axis Implants, Components and Drilling Information. BAT36d Prosthetic Flowchart. Ø6.00mm BBBT12d & BBBT24d Co-Axis Implant, Components and Drilling Information. 	Page 24 Page 25 Page 26 Page 27 Page 28 Page 29 Page 30 Page 31	
Plat	Ø4.00mm IBT12d Co-Axis Implant, Components and Drilling Information.IBT12d Prosthetic Flowchart.Ø5.00mm BAT12d & BAT24d Co-Axis Implants, Components and Drilling Information.BAT12d and BAT24d Prosthetic Flowchart.Ø5.00mm BAT36d Co-Axis Implants, Components and Drilling Information.BAT36d Prosthetic Flowchart.Ø6.00mm BBBT12d & BBBT24d Co-Axis Implant, Components and Drilling Information.BBBT12d and BBBT24d Prosthetic Flowchart.	Page 24 Page 25 Page 26 Page 27 Page 28 Page 29 Page 30 Page 31	
Plat Inst Sou Sou Pre-	Ø4.00mm IBT12d Co-Axis Implant, Components and Drilling Information. IBT12d Prosthetic Flowchart. Ø5.00mm BAT12d & BAT24d Co-Axis Implants, Components and Drilling Information. BAT12d and BAT24d Prosthetic Flowchart. Ø5.00mm BAT36d Co-Axis Implants, Components and Drilling Information. BAT36d Prosthetic Flowchart. Ø6.00mm BBBT12d & BBBT24d Co-Axis Implant, Components and Drilling Information. BBBT12d and BBBT24d Prosthetic Flowchart. form Shifting with Southern Implants Externally Hexed Implants.	Page 24 Page 25 Page 26 Page 27 Page 28 Page 29 Page 30 Page 31 Page 32	
Plat Inst Sou Sou Pre- SFI Clea	Ø4.00mm IBT12d Co-Axis Implant, Components and Drilling Information IBT12d Prosthetic Flowchart. Ø5.00mm BAT12d & BAT24d Co-Axis Implants, Components and Drilling Information BAT12d and BAT24d Prosthetic Flowchart. Ø5.00mm BAT36d Co-Axis Implants, Components and Drilling Information BAT36d Prosthetic Flowchart. Ø6.00mm BBBT12d & BBBT24d Co-Axis Implant, Components and Drilling Information BBBT12d and BBBT24d Prosthetic Flowchart. form Shifting with Southern Implants Externally Hexed Implants. rument Trays. them Screws, Alternative Impression Screws & Drill Kits. thern's Enhanced Surface. cision Attachments & Bars. Bar & Accessories.	Page 24 Page 25 Page 26 Page 27 Page 28 Page 29 Page 30 Page 31 Page 32 Page 33-35 Page 36 Page 37 Page 38 Page 39	













































External Hex Implants

Southern's range of Externally Hexed Implants are available in diameters 3.25mm, 4.0mm, 5.0mm and 6.0mm. Diameters of 7.00mm, 8.00mm and 9.00mm are also available in the External Hex MAX Implants. Implant lengths vary from 6mm to 20mm and all implants are surface roughened with Southern's extensively proven and well documented enhanced surface.


The External Hex range has an impressive prosthetic selection to accommodate a vast range of prosthetic requirements that may arise. This, together with our ideal to offer compatibility with other major brands, gives the patient more options of treatment.

The implants are made from ASTM-F67-95 Grade 4 pure titanium, with a tensile strength of 550 MPa. The surface is enhanced with abrasion and chemical conditioning. The surface has been proven by way of extensive animal and clinical trials and has been in use for more than 14 years.

IBN & IBNT Implants




Diameter 3.25mm Implants and Components




IBN Straight Site Preparation Sequence




IBNT Tapered Site Preparation Sequence






INPi Implants







INPi Drills



Site Preparation Sequence





An implant is a means to achieve a prosthetic result.

The INPi (internal drive) implant, although similar to the IBN & IBNT narrow diameter implants, adds another prosthetic option: that of Procera[™]. It has the "narrow platform" Interface and M1.6 prosthetic screws to facilitate this. In addition to this, it features Southern's proven rough surface, which has 10 years clinical documentation.



IBS / IBi, IBT / IBTi & IBPS Implants

Diameter 4.0mm Implants and Components

Cover Screw

Healing Abutments

Two-Part Healing Abutments















IBS / IBi Implants are available in: 7, 8.5, 10, 11.5, 13, 15, 18 & 20 mm. IBS Implants are pre-mounted.

IBT Implants are available in: 6, 8.5, 10, 11.5, 13, 15 & 18mm. IBT Implants are pre-mounted.

IBTi Implants are available in: 10, 11.5, 13 & 15mm.



IB Straight Site Preparation Sequence

















IBT Tapered Site Preparation Sequence






Step 1: Pilot Drill to full depth of implant

Step 2: Dedicated Drill to full depth

Step 3: Optional - Tap for cortical bone

Step 4: Place Implant

(illustration is for a 13mm implant) drill markings are at 7, 10, 13 and 15mm




BA / BAi & BAT / BATi Implants


Diameter 5.0mm Implants and Components

Two-Part Healing Abutments Cover Screw Healing Abutment SCAU5 TBA XBA **WBA** T6BA Т7ВА T5BA Ø5.5 Ø6.5 Ø7.5 Ø5.5 Ø6.5 Ø7.5 3/4/6/12 2/3/4/6 2/3/4/6/8 2/3/4/6 3/4/6 3/4/6 lengths lengths lengths **lenaths** lenaths lengths Ø5.00mm Ø5.00mm Hex 2.7mm Hex 2.7mm 0.70mm 0.70mm Enhanced Surface Enhanced BA Implants are available in: BAT Implants are available in: 6, 7, 8.5, 10, 11.5, 13, 15 & 18mm. 6, 8.5, 10, 11.5, 13, 15 & 18mm. BA Implants are pre-mounted. BA Implants are pre-mounted. Surface BAi Implants are available in: **BATi** Implants are available in: 7, 8.5, 10, 11.5, 13, 15 & 18mm. 10, 11.5, 13, 15 & 18mm. Ø5.00mm Ø4.7mm

BA Straight Site Preparation Sequence



BAT Tapered Site Preparation Sequence




BA & BAT Prosthetic Flowchart Temps Impression Laboratory **Prosthetic** Copings Components Analogues SBA 1 (Hex) GC-EX-50 (Hex) TCBA 1h (Hex) TCBA 5h (Hex) SBA 5 (Non-Hex) GC-NX-50 (Non-Hex) TCBA 1hh (Non-Hex) TCBA 5hh (Non-Hex) SCAU5 CBAU Pick-up UCLA 2 Series Screws Gold Titanium Plastic PEEK DBA DBAS12 DBAS24 CBAU-M Pick-up TITANIUM Abutments LSA12 3 Series Screws XBA (one-part) T6BA (two-part) BA & BAi Direct 20° / 35° 2/35/5 Scalloned Scalloned CERAMIC Abutments 2 Series Screws WBA (one-part) T7BA (two-part) BAT & BATI SBA16 (Hex) SBA-17-TT (Non-Hex) CIA-EX-50 (Engaging) PASSIVE Abutment SCANNING Abutment OR CBA75 Transfer Conical Abutment HBA7 TEBA1 1/2/3 CMCW1 Pick-up LSMCW1 TMCW PMCW1 GMCW1 HMCW HMCTW9 Compact Conical Abutment OR ABAMC OR 1 Series Screws 4/6 Meta**l** 4/6 Metal 1/5 Titanium CMCW2 Transfer Plastic Gold Passive Abutment тмс CMC1 Pick-up PA-MC-48 LSMC1 GMC1 нмс HMCT7 ABAMC17D-3 or ABAMC30D-4 OR OR OR Titanium 1/5 4/6 Meta**l** TMCSL Gold Passive Abutment Shouldered Abutment LT7/20 LSDBAN4 SB-DBAN1 DBAN 2/4 **Overdenture Abutment**


BBBS / BBBi & BBBT / BBBTi Implants

<u>/\</u>


Diameter 6.0mm Implants and Components




BBB Straight Site Preparation Sequence



BBBT Tapered Site Preparation Sequence



















Immediate implant placement into a molar extraction socket.

Background

The immediate placement of a conventional dental implant into a molar extraction socket poses a number of difficulties. Most significantly, the size and shape of the multi-rooted molar socket is not suited to optimal placement of a typical dental implant, often resulting in compromised implant positioning, poor primary stability or the inability to place an implant at all. This may result in the need for a waiting period of 3 to 4 months to allow for healing of the socket prior to attempting implant placement. Often, the healed extraction site presents with reduced bone height, inadequate for implant placement, resulting in the need for bone augmentation procedures, especially in the maxilla. This leads to further lengthening of treatment time with increased cost and complexity. An alternate approach has been to place a 6.00mm diameter implant into one socket of a multi-rooted extraction site, typically the palatal socket of a maxillary molar. Problems associated with the latter approach include adverse biomechanical forces resulting from the implant being off- centre and off-axis to the application of load, poor emergence profile and difficult hygiene maintenance resulting from the unavoidable buccal overhang of the restoration.

The MAX Concept

The concept of the MAX implant provides for a design of implant and a surgical protocol which makes immediate placement of the implant into a multi-rooted molar socket attainable, thus obviating the multiple problems discussed above.

The MAX implant features a body with larger than conventional diameter to fill the molar site and achieve primary stability from engagement of the perimeter bony wall of the specially prepared molar socket. The greater taper of the implant body allows for maximal preservation and engagement of inter-radicular bone within the socket of a molar with divergent roots. In the case of a molar tooth with tapering root form, the implant body has a natural fit to the socket shape. The tapered geometry of the implant allows excellent primary stability to be achieved by the threaded implant.

Surgical protocol

The ability to immediately place a MAX implant into a fresh molar extraction site represents the major advantage of this innovative treatment modality. The modality is, however, critically dependant on the preservation of the perimeter bony walls of the socket at extraction. In the case of a multi-rooted molar tooth, the extraction is best carried out by sectioning the tooth to allow removal of the roots individually, avoiding fracture of the buccal plate. If the crown of the molar is cut off horizontally, the roots can be separated and the inter-radicular bone within the socket can be removed to provide space into which the roots can be elevated. Once the roots are removed, further preparation of the socket is carried out to create a suitable tapered shape to receive the implant.

Due to the highly variable anatomy of molar roots, the operator needs to adapt the method of surgical preparation of the socket according to circumstances. Socket preparation may include all or part of the following sequence:

- a. Removal of inter-radicular bone within the socket using surgical burs or Rongeur type instruments, to create a platform for pilot hole preparation.
- b. Preparation of a pilot hole in the center of the socket or desired axis of placement of the implant. Pilot hole preparation may include the use of tapered implant drills of 4, 5 or 6 mm diameter.
- c. Shaping of the socket using a dedicated osteotome tool designed for the MAX implant which achieves deformation of the bony walls rather than further bone removal.
- d. Shaping of the socket using dedicated MAX implant drills.

The threaded MAX implant is screwed into the prepared site to achieve optimal primary stability. It is recommended that bone removed from the socket by the preceding preparation be harvested in a suction trap and utilized to fill any remaining voids around the seated implant. The occurrence of voids is frequent due to the irregular shape of molar sockets.

Restorative Protocol

The restorative interface of the MAX implant features a wide platform and external hex.

The implant is restored by "platform switching" i.e. the platform of the restorative component has a smaller diameter than that of the implant platform. Restorative procedures are otherwise similar to those of other externally hexed implants. See flow charts of available restorative options.

Advantages of the MAX protocol

- 1. Immediate placement of the MAX implant into a molar socket has the following advantages:
 - -preservation of alveolar bone
 - -avoidance of separate grafting procedures
 - -shortened treatment time for molar replacement
 - -reduced number of surgical procedures
 - -reduced cost and complexity of treatment
- 2. Excellent primary stability is achieved by the tapered, threaded design engaging perimeter bony walls of the socket.


















MAX-7 Implants




Diameter 7.0mm Implants used with diameter 5.0mmComponents

drill markings are at 7, 10, 13 and 15mm

Two-Part Healing Abutments Cover Screw Healing Abutments SCAU5 TBA XBA WBA T6BA T7BA T5BA Ø5.5 Ø6.5 Ø7.5 Ø5.5 Ø6.5 Ø7.5 2/3/4/6/8 3/4/6/12 2/3/4/6 2/3/4/6 3/4/6 3/4/6 lengths lengths lengths lengths lengths lengths Ø5.5mm Hex 2.7mm 0.70mm Ø6.0mm MAX-7 are pre-mounted and are Surface Enhanced available in the following lengths: 7mm code: MAX-7-7 9mm code: MAX-7-9 Thread Pitch code: MAX-7-11 11mm 0.8mm Ø7.00





D-70TP-11 D-TAP-MAX7-11

D-30T-M15

D-60TP-11

D-RB-MS

D-20T-M15


MAX-7 Prosthetic Flowchart Temps Impression Laboratory **Prosthetic** Copings Components Analogues SBA 1 (Hex) GC-EX-50 (Hex) TCBA 1h (Hex) TCBA 5h (Hex) SBA 5 (Non-Hex) GC-NX-50 (Non-Hex) TCBA 1hh (Non-Hex) TCBA 5hh (Non-Hex) CBAU Pick-up UCLA 2 Series Screws SCAU5 Gold Titanium Plastic PEEK DBAS12 DBA DBAS24 TITANIUM Abutments CBAU-M Pick-up LSA12 3 Series Screws XBA (one-part) T6BA (two-part) Direct 20° / 35° 2/3.5/5 Scalloned CERAMIC Abutments 2 Series Screws WBA (one-part) T7BA (two-part) SBA16 (Hex) SBA-17-TT (Non-Hex) CIA-EX-50 (Engaging) PASSIVE Abutment SCANNING Abutment OR CBA75 Transfer Conical Abutment HBA7 TEBA1 1/2/3 CMCW1 Pick-up LSMCW1 TMCW PMCW1 GMCW1 HMCW HMCTW9 Compact Conical Abutment OR ABAMC OR 1 Series Screws 4/6 Meta**l** 4/6 Metal 1/5 Titanium CMCW2 Transfer Plastic Gold Passive Abutment тмс CMC1 Pick-up PA-MC-48 GMC1 нмс HMCT7 ABAMC17D-3 or ABAMC30D-4 OR OR OR Titanium 1/5 4/6 Meta**l** TMCSL Gold Passive Abutment Shouldered Abutment LT7/20 LSDBAN4 SB-DBAN1 DBAN 2/4 **Overdenture Abutment**

MAX-8 Implants



Diameter 8.0mm Implants used with diameter 6.0mmComponents






MAX-8 are pre-mounted and are available in the following lengths:


7mm code: MAX-8-7 9mm code: MAX-8-9 11mm code: MAX-8-11



MAX-8 Site Preparation Sequence



Earlier revisions of the MAX drills are 2.4mm longer than the implant. These can easily be identified by the lazer marking on the body of the drill. Current drills are marked on the shank.



MAX-9 Implants



Diameter 9.0mm Implants used with diameter 7.0mmComponents



Healing Abutments

TB9MAX


Two-Part Healing Abutment












MAX-9 are pre-mounted and available in the following lengths:

7mm code: MAX-9-7 9mm code: MAX-9-9 11mm code: MAX-9-11


MAX-9 Drills & Additional Instrumentation







D-TAP-MAX9-7 / 9 / 11




MAX-9-PG-7 / 9 / 11

Profile Gauges





MAX-9 Site Preparation Sequence



Step 1: Pilot Drill to full depth of implant

Step 2: Final Twist Drill to full depth

Step 3: Dedicated Drill to full depth


Step 4: Dedicated Drill to full depth

Step 5: Dedicated Tap for cortical bone

Step 6: Place Implant

tilmm (illustration is for a 11mm implant) drill markings are at 7, 10, 13 and 15mm

Earlier revisions of the MAX drills are 2.4mm longer than the implant. These can easily be identified by the lazer marking on the body of the drill. Current drills are marked on the shank.



MAX Ø10.00mm Implants are also available in lengths 7, 9, 11 and 13mm as "special items". For more information please contact Southern Implants Head Office.

















The Co-Axis implant is indicated for use in situations where the long axis of a conventional implant would not coincide with the long axis of the restoration and would therefore result in a restorative compromise.

The most common example of this is encountered where an implant is placed in the anterior maxilla at a labially inclined angle, as dictated by the anatomy of the alveolus, resulting in the screw access hole of the prosthetic crown passing through the labial face of the crown. The Co-Axis implant effectively solves this problem by having the prosthetic platform and screw hole of the implant tilted at an angle of 12 or 24 degrees to the long axis of the implant. The axis of the retaining screw is therefore also offset within the body of the implant.

The Co-Axis concept can be applied to solve many other situations where inclined placement of implants is either unavoidable or even an advantage. For example where avoidance of anatomical structures dictates (eg maxillary sinus, mental foramen) or where bony anatomy can be maximised by inclined placement of an implant. An elegant and truly innovative solution to a frequent problem in implant dentistry.

Advantages

- 1. The Co-Axis solution greatly simplifies the restorative treatment of an inclined implant by eliminating the need for angle correcting abutments or custom abutments. This reduces the number and cost of components required, reduces the complexity and cost of laboratory work as well as the number of patient visits required.
- 2. Esthetic advantages result from having no need for labially placed screw access holes.
- 3. Avoidance of anatomical structures by inclined implant placement, without incurring prosthetic complications, is made possible by exploiting the Co-Axis concept.
- 4. The Co-Axis implant allows for maximal utilisation of available bone anatomy and may result in the advantage of being able to use either a longer or larger diameter implant.
- 5. Screw retained restorations can be used instead of cemented restorations
- 6. Immediate loading protocols are greatly facilitated by use of the Co-Axis implant by making screw retained restorations routinely attainable.
- 7. The tapered Co-Axis implant provides an anatomically correct implant for ideal use in the anterior maxilla.

Placement

Implant site drilling is done in the conventional manner, using tapered drills, but in such a way that the 12 or 24 degree angle correction that the implant head offers, is exploited. This allows for best utilisation of available bone. The bone preparation is best done to a depth slightly greater than the intended implant length e.g. drill to 15mm for a 13mm implant. This technique will assure full seating of the implant to the desired depth.

The use of a Co-Axis direction indicator in the prepared site after pilot drilling and after final drilling, will allow the operator to visualise the restorative axis of the implant, before committing to implant placement. The implant is mounted on a special fixture mount which compensates for the tilted implant head and thus results in straight line placement of the implant, just as any other fixture mounted implant.

The fixture mount is marked with a clearly visible dimple which is used as a reference of the rotational position of the implant relative to the angle correction ie. the platform of the implant is angled AWAY FROM the dimple. This means that, in a typical anterior maxillary placement, the dimple should end up being positioned on the mid labial aspect of the site. In the event that the fixture mount and implant need to be reassembled, care must be taken to align the fixture mount dimple with the notch cut into the top of the implant hexagon. This will ensure concentric rotation on the assembly.

The implant should be placed to a depth dictated by a number of factors, including esthetic demand for concealment of the implant collar, thickness of soft tissue and personal preferences of the operator or restorative dentist. Typically, if the broad collar on the labial aspect of the implant is placed to the crest of the labial bone height, it may result in other areas of the collar being sub-crestal. This would be an indication for using a bone mill before seating the preferred abutment.

The Co-Axis implant is specifically designed with a finer thread pitch than most other implants, with the result that a full rotation of the implant results in a depth change of only 0.6mm. This feature makes it easier to achieve the correct positioning of the 12 degree angle (dimple indicator) in combination with the correct depth of placement.


















IBT12d Co-Axis Implants







IBT12d are pre-mounted and are available in the following lengths:

 8.5mm
 code:
 IBT12d-8.5

 10mm
 code:
 IBT12d-10


 11.5mm
 code:
 IBT12d-11.5

 13mm
 code:
 IBT12d-13


 15mm
 code:
 IBT12d-15


 18mm
 code:
 IBT12d-18

IBT12d Additional Instrumentation




IBT12d Site Preparation Sequence






BAT12d & BAT24d Co-Axis Implants



BAT12d / 24d Additional Instrumentation



BAT12d / 24d Site Preparation Sequence



BAT12d & BAT24d Prosthetic Flowchart Prosthetic Temps Impression Laboratory Copings Components Analogues SBA 1 (Hex) GC-EX-50 (Hex) TCBA 1h (Hex) SBA 5 (Non-Hex) GC-NX-50 (Non-Hex) TCBA 1hh (Non-Hex) SCALIS CBAU Pick-up UCLA 2 Series Z Screws Gold Titanium PEEK DBA TITANIUM Abutment CBAU-M LSA12 3 Series Z Screws XBA (one-part) T6BA (two-part) BAT12d Direct 2/3.5/5 CER-ZR-58 CERAMIC 2 Series Z Screws WBA (one-part) T7BA (two-part) BAT24d SBA16 (Hex) SBA-17-TT (Non-Hex) CIA-EX-50 (Engaging) PASSIVE Abutment SCANNING Abutment CBA75 Transfer Uses 3 Series Screws PA-MC-60 CMCW1 Pick-up **Compact Conical Abutme** LSMCW1 PMCW1 GMCW1 TMCW HMCW OR 1 Series Screws 4/6 Metal 4/6 Metal 1/5 Titanium CMCW2 Transfer Plastic 1/2/3/4/5 Passive Abutment Shouldered Abutment SB-DBAN1 Waxing sleeve for uncut DBAN LSDBAN4 P-CAP-W40 2/4 Analogue for uncut DBAN Plastic

NOTE: 1. Platform Shifting is possible and in some cases the preferred restorative protocol for this implant. For Platform Shifting use the restorative options of the IBT12d (page 25).

2.The BAT12d range can use the same prosthetics as for the BA & BAT Implant range. (See page 11)





This innovative concept was specifically designed with the tilted Implant techniques in mind.

By making the angle of the implant head 24d, complicated prosthetic procedures where angled abutments are needed are eliminated. The restorative platform of this implant is $\emptyset 5.00$, the same as the BA range of implants. Special prosthetics are required, as depicted above, due to the internal thread of the implant being slightly shallower. The Z Series screws must be used.

BAT36d Co-Axis Implants



Diameter 5.0mm Implants used with Diameter 4.0mm Components


Cover Screw

Healing Abutments











BAT36d are pre-mounted and are available in the following lengths:

code: BAT36d-10 10mm 11.5mm code: BAT36d-11.5 13mm code: BAT36d-13 15mm code: BAT36d-15 18mm code: BAT36d-18

BAT36d Additional Instrumentation







D-TAP-L-50


Please note:

Point 1

This corner of the drill is to be at bone level.

Point 2

This corner of the drill will be subcrestal.



D-52TP-13

D-50TP-13 (soft bone)

Step 1: Pilot Drill to full depth of implant Step 2: Dedicated Drill to full depth

This is a Ø5.1 implant and hence the gold coloured drill D-52TP-13 is the

ideal final drill. However, in soft bone, D-50TP can be used to obtain greater primary stability.

Step 3: Optional - Tap for cortical bone In hard bone, tapping of the cortical bone is recommended.


Step 4: Place Implant & remove fixture mount

(illustration is for a 13mm implant)

drill markings are at 7, 10, 13 and 15mm

D-20T-M15


D-RB-MS



BBBT12d & BBBT24d Co-Axis Implants












BBBT12d / 24d Additional Instrumentation



BBBT12d / 24d Site Preparation Sequence





Platform Shifting

The platform shifting concept implies the use of prosthetic components having a platform diameter undersized when compared to the diameter of the implant platform. In this way, the prosthetic connection is displaced horizontally inwards from the perimeter of the implant platform, creating an angle or step between the abutment and implant. This concept has been widely published with claims of better bone levels.

Important: For platform shifting to be effective, it is important to carry out the under sizing of the components during all phases of the implant treatment, i.e. from placement of the implant through to final restoration.

The principles and requirements for improved hard and soft tissue response by way of platform shifting are now well established in the literature. Of importance is the establishment of the shifted dimension from the time of implant placement, through all restorative phases.

Platform Shifting with Ø4.0 Externally Hexed Implants



The Ø4.0 Implants get restored with Standard Platform Prosthetics as per page 7 to apply platform shifting.




Platform Shifting with Ø5.0 Externally Hexed Implants

BA **BAT** Ø5.0 BAT12d BAT24d



The Ø5.0 Implants get restored with Standard Platform Prosthetics as per page 7 for the BA, BAT and BAT12d implants. The BAT24d gets restored using the Co-Axis Standard Platform Prosthetics as per page 25 to apply platform shifting.



Platform Shifting with Ø6.0 Externally Hexed Implants


BBBS $\emptyset 6.0$ **BBBT** BBBT12d BBBT24d



The Ø6.0 Implants get restored with Ø5.0 Prosthetics as per page 9 to apply platform shifting. The BBBT24d gets restored using the Co-Axis Ø5.0mm Platform Prosthetics as per page 29 to apply platform shifting.




I-HEX-U






I-HEX-1




I-DRILL




These trays fit into the I-CASSETTE-2 and are autoclavable. Please see CAT-1039 for sterilizing and cleaning instructions.

I-MAX




I-PROS-B



These trays fit into the I-CASSETTE-2 and are autoclavable. Please see CAT-1039 for sterilizing and cleaning instructions.

Southern Implants Screws



Series 1 Screws (M1.4)

10-15Ncm Head Diameter 2.25mm

Series 2Z & 2 Screws (M2)

32-40Ncm

Head Diameter 2.70mm

* = Long Versions Available

Series 3Z & 3 Screws (M2)

32-40Ncm


Head Diameter 2.45mm

* = Long Versions Available



Please Note: Series 2Z Screws have two grooves on the head for easy recognition.

Alternative Impression Screws



Drill Kits

















Southern's Enhanced Surface


The Southern enhanced surface is not a "coating", it is an **abraided rough surface of Rutile Titanium**. This is the same dense form of titanium common to "machined" surface implants. (The anodic oxidation surfaces are not Rutile Titanium; they are a mixture of anatase and amorphous titanium which are less dense and softer forms of titanium).


- A. The first experimentation with this Southern Enhanced surface was in **1992**. After extensive validation it was **put into widespread clinical use in 1997**. It is achieved by a subtractive process in which specifically sized and shaped, sharp cornered, Alumina particles (Al₂0₃) are blasted with decontaminated air onto the implant surface with controlled pressure, displacement and time. Every batch of Al₂0₃ particles are subject to SEM analysis to ensure consistent shape and size.
- B. The particle size of 110μm is supported by the work of Soskalne (Israel) and Wennerberg (Sweden) on the one hand and Ronald (Norway) on the other. Based on their research, greatest bone to titanium bond strength is obtained with abrasion particles **greater that 75μm and less than 170μm**. This partly explains why the 250μm particles used to produce the SLA surface, need a secondary acid etching process as the 250μm particles alone, leave too course a surface. Astra wisely use 110μm particles for producing their rough surfaces, but their particles are of titanium oxide and not Alumina.
- C. Szmukler-Moncler has analyzed and compared the popular implant surfaces in publications and a presentation at the AO, San Francisco 2004. He reports that the Southern Surface is **remarkably consistent** and **free of containments** whilst those that are acid etched or oxidized are shown to be highly variable. It is extremely difficult to control acid etching and oxidation in an industrial manufacturing process. This is one reason why Southern does not use acid etching or anodic oxidation.
- D. There seems to be consensus in the literature that "moderately rough" surfaces have no great risks for the patient and are therefore **safe to use**. Moderately rough was defined by Albrektsson as S_a = 1.0 to 2.0 μ m (applied Osseointegration Research Vol 5, 2006) and our surface has S_a = 1.43 in one published study and S_a = 1.55 on implants recently analyzed by Prof Ann Wennerberg in 2006.

Dr Mats Wikström, Chief of Clinics, Branemark Center Goteborg, in 2007 concluded that the Southern surface is **one of the three best documented moderately rough surfaces on the market**.


Prof Alan Payne, Oral Implantology Research Group, University of Otago, is conducting Randomized Clinical Trials (RCTs) involving Southern Implants' rough surface. 2008 signified the 10 year follow-up in the mandibles and 7 year follow-up in the Maxilla. The 8 year and 5 year results are published in Cochrone Collaboration reports. Standardized radiographs show marginal bone loss of all the implants to be well within the criteria set by Albrektsson & Zarb 1993, 1998 as well as Fourmousis & Bragger 1999.








In conclusion, it is a well documented surface with a consistent manufacturing process and holds extremely low risk.





SFI Bar

The SFI-Bar is the innovative bar solution for removable dentures on implants in both upper and lower jaws. The bar ensures stress-free hold of the prosthesis on the implants, thus increasing the patient's comfort. All components are prefabricated. Therefore problems associated with soldering and welding conventional bars are eliminated completely.


Moreover due to it's simplicity the SFI-Bar can be easily and efficiently adapted and placed aspiration-safe to the individual patient's mouth situation chair side. Its functional design ensures great flexibility and is indicated in most cases. With the appropriate abutments, the SFI-Bar can be used with almost any implant system.

The SFI-Bar 2-implant and 4-implant are the standard versions. These can be upgraded with the SFI-Bar Add-on Kit (0500 0668) to solutions for 3, 5 and 6 implants.



- * Tension free, excellent and stable fit.
- * Simply ingenious, thanks to telescope-like connections and individual shortening.
- * Possible to fit the SFI-Bar 2-implant directly in the mouth.





The SFI Bar System and all it's accessories are only sold through Southern Implants in South Africa. For more information call Southern Implants.

Cleaning & Sterilization Procedure Guidelines

Components	Surgical Drills	Surgical and Prosthetic Tools	Surgical and Prosthetic Trays	Torque Wrenches / Ratchet			
Components							
Warnings	Never use blunt or damaged tools.	Never use blunt or damaged tools.	Do not expose to temperatures higher than 140°C.	Do not expose to temperatures higher than 160°C.			
Never store instruments while they are still wet or moist. This will result in corrosion and degradation of cutting edges. Inspect all instruments							
Limitations on Reprocessing	Taper Drills and Pilot Drills - 40 times. Trephine Drills - 10 times.	Drivers and General Instruments - Repeated processing has minimal effect. Bone Mills - 15 times. Taps, Tissue Punch, Osteotomes, counter bores, Thread cutters - 40 times.	Repeated processing has minimal effect on these trays.	500 times or 3 years, whichever is reached first (service can extend life, but not indefinately).			
	Instructions						
Point of use:	Remove excess soil with running water.	Remove excess soil with running water.	Remove excess soil with disposable cloth / paper wipe.	Remove excess soil with running water.			
Containment and Transportation	No Particular Requirements: It is recommended that instruments are cleaned and dried as soon as it is reasonably practical following use.						
Preparation for cleaning	Disassembly from Handpiece.	Disassembly from Handpiece.	No Particular Requirement.	Remove any connecting parts.			
Cleaning: Automated	Rinse with luke-warm water for 3 minutes. Remove the hardened soil with a soft bristle brush. Sonicate for 20 minutes in an ultrasonic cleaner using a detergent suitable for surgical instruments (SteriTech Instrument Cleaner - 5% dilution). Rinse in running water and dry.	Rinse with luke-warm water for 3 minutes. Remove the hardened soil with a soft bristle brush. Sonicate for 20 minutes in an ultrasonic cleaner using a detergent suitable for surgical instruments (SteriTech Instrument Cleaner - 5% dilution). Rinse in running water and dry.	Rinse in luke-warm water and remove soil with a soft bristle brush and enzymatic detergent. Rinse off with luke-warm water and dry.	Rinse with luke-warm water for 3 minutes. Remove the hardened soil with a soft bristle brush. Sonicate for 20 minutes in an ultrasonic cleaner using a detergent suitable for surgical instruments (SteriTech Instrument Cleaner - 5% dilution). Rinse in running water and dry.			
Cleaning: Manual	Manual cleaning is not practical and is therefore not advised.						
Disinfection	If detergent is not a high level disinfectant: Sonicate / rinse with 70 % ethanol for 5 minutes.	If detergent is not a high level disinfectant: Sonicate / rinse with 70 % ethanol for 5 minutes.	Trays can be wiped / sprayed down with 70 % ethanol.	If detergent is not a high level disinfectant: Sonicate / rinse with 70 % ethanol for 5 minutes.			
Drying	Allow components to dry completely before sterilizing.						
Maintenance	Damaged or blunt instruments should not be used.			Add one drop of handpiece oil to the internal mechanism after autoclaving.			
Inspection and further Testing	Inspect all instruments visually for damage and bluntness.			Check whether the ratchet mechanism is working smoothly.			
Packaging	Single: A standard packaging material may be used. Ensure that the pack is large enough to contain the instrument without stressing the seal. (Recommendation: Use a packaging system that conforms to ISO 11607).						
	In sets: Instruments can be loaded into dedicated instrument trays before sterilization. Place the tray in an appropriate sterilization bag.						
Sterilization	Autoclave at 121 ° C (250 ° F) for a minimum of 30 minutes.						
Storage	Never store instruments while they are still wet or moist. This will result in corrosion and degradation of cutting edges.						
Additional Information When sterilizing multiple instruments in one autoclave cycle, ensure the the sterilizer's maximum load is not exceeded. Refer to sterilizer's instructions for use. Ascertain the size of the instrument trays prior to sterilization.							
	<u> </u>						







Complimentary Manuals & Instructions:

004
005
022
023
032
025
006
001
036
017
018
800
083
010
024
038

Labeling Symbols:

The following symbols are used on our packaging labels and they indicate the following:

1: "Use by"

LOT

2: "Batch code"

(2)

3: "Do not reuse"

STERILE R

4: "Sterilization using Irradiation"

Δ

5: "Caution"

6: "Consult instruction for use"

i

7: CE mark







www.southernimplants.com

Head Office: South Africa Southern Implants (Pty) Ltd. Tel: +27 12 667 1046 Fax: +27 12 667 1029 info@southernimplants.com Head Office: Americas / Asia Southern Implants Inc. Tel: +1 949 273 8505 Fax: +1 949 273 8508 info@southernimplants.us

United Kingdom

Southern Implants UK Tel: +44 208 998 0063 Fax: +44 208 997 0580 hilary@southern-ventura.co.uk

S.I.R srl Tel: +39 04 556 8351 Fax: +39 04 556 8351 e.caprini@southernimplants.it

Australia Henry Schein I Halas Dental Tel: +61 2 9697 6288 Fax: +61 2 9697 6250 info@henryschein.com.au

Spain / Portugal

Contactodent Tel: +351 214 693 332 Fax: +351 214 693 329 southernimplants@sapo.pt

New Zealand Southern Implants Ltd NZ Tel: 0800 246 752 Fax: +64 9 430 2836 dkshep@xtra.co.nz

Nordic Countries

Protera AB Tel: +46 31 291078 Fax: +46 706 150078 villy@protera.se

Southern Implants
Tel: +49 7121 490 620
Fax: +49 7121 491 717
info@southernimplants.de

Turkev

Ekodent Tel: +212 343 5233 selvasoyturk@ekodent.com www.ekodent.com

Benelux ProScan bvba Tel: +32 11 822 650 Fax: +32 11 822 651 info@proscan.be

Namibia

Skydancer Tel: +64 61 225 152 Fax: +64 61 235 630 mfdental@iway.na

Southern Implants Tel: +30 210 898 2881 Fax: +30 210 595 2543 info@southernimplants.gr

Russia

Implant Angle Ltd Tel: +7 495 649 1732 Fax: +7 495 789 9133 info@implantangle.ru